



Abbott Analytical



Consulting Scientists to the Disinfectant Industry

Certificate of Analysis

Sample(s): One sample of Remuvit Concentrate

Received from: Klenitise Ltd. 14 Cherry Grove, Sketty, Swansea, SA2 8AS

Date received: 16 May 2011 **Date tested:** 19 May 2011

Certificate no: 11E.040Sh-KR.CLE **Certificate date:** 23 May 2011

Sample ref: 11E/040 **Page:** 1 of 2

Analysis required: EN 1650, Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of fungicidal or yeasticidal activity of chemical disinfectants and antiseptics used in food, industrial, domestic and institutional areas - Test method and requirements (phase 2, step 1)

Product stored at: Room temperature

Active substance: Not declared

Test conditions: Dirty

Interfering substance: 3.0g/l bovine albumin

Product test concentration: 20% v/v

Product diluent used during test: Sterile hard water 300mg/l CaCO₃

Contact time: 15 minutes

Test temperature: 20°C ± 0.5°C

Neutralising solution: 30g/l polysorbate 80, 3g/l lecithin, 1g/l histidine, 1g/l cysteine

Incubation temperature: 30°C ± 1°C

Identification of fungal/yeast strain(s) used: *Stachybotrys chartarum* NCPF 7587

D C Watson



Abbott Analytical



Consulting Scientists to the Disinfectant Industry

23 May 2011

Certificate No: 11E.040Sh-KR.CLE

Page: 2 of 2

Test results:

Test Organism	<i>Stachybotrys chartarum</i>	
Validation Suspension (N_v)	Vc1 326	Vc2 304
	$\bar{x} = 315$	
Experimental Control (A)	Vc1 288	Vc2 312
	$\bar{x} = 300 \geq 0.5N_{v_0}$	
Neutraliser Control (B)	Vc1 264	Vc2 296
	$\bar{x} = 280 \geq 0.5N_{v_0}$	
Method Validation (C)	Vc1 274	Vc2 310
	$\bar{x} = 292 \geq 0.5N_{v_0}$	
Test Suspension	10^{-5} Vc1 320	Vc2 308
	10^{-6} Vc1 38	Vc2 31
(N)	$\bar{w} = 3.17 \times 10^7$	
	lg N = 7.50	
(N₀ = 0.1N)	lg N ₀ = 6.50	
Results	10^0 Vc1 <14	Vc2 <14
	(Na) $10\bar{x} < 140$	
	lg Na < 2.15	
(R)	lg R > 4.35	
Pass: lg R \geq 4	PASS	

Vc = plate count per ml
 \bar{x} = average of Vc1 and Vc2
 \bar{w} = weighted mean of \bar{x}
R = reduction (lg R = lg N₀ - lg Na)

Conclusion:

This batch of Remuvit Concentrate, when diluted to 20% v/v, passes the requirements of EN 1650 for fungicidal/yeastocidal activity in 15 minutes at 20°C under dirty conditions against the reference organism detailed.

D C Watson

D C Watson BSc, CBiol, MSB, MIFST, ACIEHO
PO Box 95, New Ferry, Wirral, CH62 6HA
Tel: 0151 637 3331 Mob: 07767 871275
email: abbottanalytical@hotmail.co.uk